

Product Specification

- 2000ppm hypochlorous acid in purified water
- Ready-to-use with 2 year shelf life
- 0.2 micron filtered under Grade A air in a Grade B/C cleanroom
- Double bagged in polyethylene linear tear bags
- Trigger spray with jet or spray option
- Relevant certificates of analysis and sterility provided with every batch.
- Full supporting documentation information is available on request.

Efficacy

Test	Description	Log Reduction	Time	Test	Description	Log Reduction	Time
EN1276	<i>E. hirae</i>	>log 7	1 min	EN13704	<i>C. difficile</i>	>log 6	1 min
EN1276	<i>S. aureus</i>	>log 7	1 min	EN13704	<i>B. subtilis</i>	>log 6	1 min
EN1276	<i>P. aeruginosa</i>	>log 7	1 min	EN13704	<i>B. pumilis</i>	>3.88	1 min
EN1276	<i>E. coli</i>	>log 7	1 min	EN13704	<i>B. cereus</i>	>3.24	1 min
EN14776	Poliovirus	4.33	30 secs	EN13697	<i>B. subtilis</i>	> log 6	1 min
EN1650	<i>A brasiliensis</i>	>log 4	1 min	EN16615	<i>A brasiliensis</i>	>6.13 test area 1	1 min
EN1650	<i>C. albicans</i>	> log 4	1 min	EN16615	<i>B. subtilis</i>	>7.11 test area 1	1 min

Product Application

Ready-to-use so can be applied directly to the surface ensuring even coverage. Use an appropriate cleanroom wipe or mop and a recommended wiping technique for optimum contamination control. Leave for required contact time and wipe to dry. Contec ProChlor is low residue but where residues are an issue they are easily removed with either water or alcohol.

Part No.	Description	Size	Packaging
SBT102PC	Contec Sterile ProChlor Stabilised HypoChlorous acid in purified water, trigger spray	1L	6 bottles x case
SBC102PC	Contec Sterile ProChlor Stabilised HypoChlorous acid in purified water, capped	1L	6 bottles x case
SBC502PC	Contec Sterile ProChlor Stabilised HypoChlorous acid in purified water, capped	5L	2 bottles x case
FBT102PC	Contec ProChlor Stabilised HypoChlorous acid in purified water, trigger spray	1L	6 bottles x case
FBC102PC	Contec ProChlor Stabilised HypoChlorous acid in purified water, capped	1L	6 bottles x case
FBC502PC	Contec ProChlor Stabilised HypoChlorous acid in purified water, capped	5L	2 bottles x case



Use biocides safely. Always read the label and product information before use.